

Product Datasheet

Troponin I (41-49) Antibody HM412

Product Name	Troponin I (41-49) Antibody
Catalog Number	HM412
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified human Troponin complex
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Troponin I (epitope 41-49 a.a.)
Matched Pairs	HM413 or HM414 (capture)
Preservatives	0.09% Sodium Azide