

Product Datasheet

Troponin I (83-95) Antibody HM414

Product Name	Troponin I (83-95) Antibody
Catalog Number	HM414
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Purified human Troponin complex
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Troponin I (epitope 83-95 a.a.)
Matched Pairs	HM412 or HM413 (detection)
Preservatives	0.09% Sodium Azide