

Product Datasheet

Hepatitis B Core Antigen (HBcAg) Antibody HM415

Product Name	Hepatitis B Core Antigen (HBcAg) Antibody
Catalog Number	HM415
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Recombinant HBcAg
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	HBcAg
Matched Pairs	HM416 (capture)
Preservatives	0.05% Proclin 300