

# Product Datasheet

## Human Granzyme B ELISA Kit 850.790.048

Product Name	Human Granzyme B ELISA Kit
Catalog Number	850.790.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	20 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with Trail, TNFα, IFNγ, Fas, IL-1α, IL-4, IL-6, IL-8, IL-13
Specificity	Recognizes both natural and recombinant human Granzyme B
Sample Types	Serum, Plasma, Cell culture supernatant