

Product Datasheet

Influenza B Virus Antibody HM427

Product Name	Influenza B Virus Antibody
Catalog Number	HM427
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Immunogen	Inactivated Influenza B Virus
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Adenovirus, and RSV
Specificity	Influenza B Virus Nucleoprotein
Matched Pairs	HM428 (detection)
Preservatives	0.05% Sodium Azide