

Product Datasheet

Adenovirus Antibody HM429

Product Name	Adenovirus Antibody
Catalog Number	HM429
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Immunogen	Inactivated Adenovirus
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥ 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Influenza B, and RSV
Specificity	Hexon Protein
Matched Pairs	HM430 (detection)
Preservatives	0.05% Sodium Azide