

Product Datasheet

Human CD213a2 ELISA Kit 850.800

| | |
|---------------------|--|
| Product Name | Human CD213a2 ELISA Kit |
| Catalog Number | 850.800 |
| Description | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 156-5000 pg/ml |
| Sensitivity | 150 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | IL-13Ra2 |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity has been observed with TNF-a, TRAIL, IFN-g, Fas, IL-1a, IL-4, IL-6, IL-8, IL-13, Granzyme B |
| Specificity | Recognizes both natural and recombinant IL-13Rα2 |
| Sample Types | Serum, Plasma, Cell culture supernatant |