

Product Datasheet

Respiratory Syncytial Virus (RSV) Antibody HM437

Product Name	Respiratory Syncytial Virus (RSV) Antibody
Catalog Number	HM437
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF
Immunogen	Purified native RSV
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2a
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Influenza B, and Adenoviruses
Specificity	Fusion protein of RSV strains A & B
Matched Pairs	HM438 or HM439 (detection)
Preservatives	0.05% Sodium Azide