

# Product Datasheet

## Human CD213a2 ELISA Kit 850.800.192

Product Name	Human CD213a2 ELISA Kit
Catalog Number	850.800.192
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	156-5000 pg/ml
Sensitivity	150 pg/ml
Tested Applications	ELISA
Alternative Names	IL-13Ra2
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity has been observed with TNF-a, TRAIL, IFN-g, Fas, IL-1a, IL-4, IL-6, IL-8, IL-13, Granzyme B
Specificity	Recognizes both natural and recombinant IL-13Rα2
Sample Types	Serum, Plasma, Cell culture supernatant