

Product Datasheet

Troponin I (41-60) Antibody HM456

Product Name	Troponin I (41-60) Antibody
Catalog Number	HM456
Description	Protein A Purified and Filtered Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, TIA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified and Filtered Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	≥ 95 %
Buffer	20 mM Sodium Phosphate, pH 7.2 with 150 mM NaCl
Specificity	Troponin I (epitope 41-60 a.a.)
Matched Pairs	HM233 (capture)