

Product Datasheet

MDPV Antibody HM457

Product Name	MDPV Antibody
Catalog Number	HM457
Description	DEAE Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	MDPV-BSA
Alternative Names	Methylenedioxyprovalerone
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE and HPLC
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	MDPV
Matched Pairs	LA011 (capture)
Preservatives	0.05% Sodium Azide