

Product Datasheet

Human IL-12p70 HS ELISA Kit 850.810.048

Product Name	Human IL-12p70 HS ELISA Kit
Catalog Number	850.810.048
Description	Diaclone High Sensitivity ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.78-25 pg/ml
Sensitivity	0.75 pg/ml
Tested Applications	ELISA
Alternative Names	IL-12
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-12. Specific to p70 heterodimer.
Sample Types	Serum, Plasma, Cell culture supernatant