

Product Datasheet

Japanese Encephalitis Virus (JEV-NS1) Antibody HM467

Product Name	Japanese Encephalitis Virus (JEV-NS1) Antibody
Catalog Number	HM467
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Japanese Encephalitis Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with other Flaviviruses
Specificity	Japanese Encephalitis Virus NS1 protein
Preservatives	0.05% Sodium Azide