

Product Datasheet

Tick-borne Encephalitis Virus (TBEV-NS1) Antibody HM478

Product Name	Tick-borne Encephalitis Virus (TBEV-NS1) Antibody
Catalog Number	HM478
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Tick-borne Encephalitis Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with other Flaviviruses
Specificity	Tick-borne Encephalitis Virus NS1 protein
Preservatives	0.05% Sodium Azide