

Product Datasheet

Yellow Fever Virus (YFV-NS1) Antibody HM492

| | |
|----------------------------|--|
| Product Name | Yellow Fever Virus (YFV-NS1) Antibody |
| Catalog Number | HM492 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF, IFA |
| Immunogen | Purified Yellow Fever Virus |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1, kappa |
| Clonality | Monoclonal |
| Purity | > 90 % by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride |
| Cross Reactivity | No reactivity with other Flaviviruses |
| Specificity | Yellow Fever Virus NS1 protein |
| Preservatives | 0.05% Sodium Azide |