

Product Datasheet

Yellow Fever Virus (YFV-NS1) Antibody HM494

Product Name	Yellow Fever Virus (YFV-NS1) Antibody
Catalog Number	HM494
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, IFA
Immunogen	Purified Yellow Fever Virus
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with other Flaviviruses
Specificity	Yellow Fever Virus NS1 protein
Matched Pairs	HM495 (detection)
Preservatives	0.05% Sodium Azide