

# Product Datasheet

## Human TRAIL R2 ELISA Kit 850.830.096

Product Name	Human TRAIL R2 ELISA Kit
Catalog Number	850.830.096
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	31.25-1000 pg/ml
Sensitivity	6 pg/ml
Tested Applications	ELISA
Alternative Names	DR5, CD262
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity was observed with TRAIL, CD117, IL-6R, IL-2R, CD116, TRAIL R1, TRAIL R3, TRAIL R4, CD178 and Granzyme B.
Specificity	Specific for human TRAIL R2 / DR5
Sample Types	Serum, Cell lysates, Cell culture supernatant