

## Product Datasheet

### Human TRAIL R2 ELISA Kit 850.830.192

<b>Product Name</b>	Human TRAIL R2 ELISA Kit
<b>Catalog Number</b>	850.830.192
<b>Description</b>	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
<b>Range</b>	31.25-1000 pg/ml
<b>Sensitivity</b>	6 pg/ml
<b>Tested Applications</b>	ELISA
<b>Alternative Names</b>	DR5, CD262
<b>Brand</b>	Diaclone
<b>Storage</b>	Store at +4°C
<b>Cross Reactivity</b>	No cross reactivity was observed with TRAIL, CD117, IL-6R, IL-2R, CD116, TRAIL R1, TRAIL R3, TRAIL R4, CD178 and Granzyme B.
<b>Specificity</b>	Specific for human TRAIL R2 / DR5
<b>Sample Types</b>	Serum, Cell lysates, Cell culture supernatant