

Product Datasheet

Methicillin-resistant Staphylococcus aureus (MRSA) Antibody HM504

Product Name	Methicillin-resistant Staphylococcus aureus (MRSA) Antibody
Catalog Number	HM504
Species/Host	Mouse
Tested Applications	ELISA, IP, WB, IHC
Immunogen	Purified MRSA protein (mecA)
Brand	EastCoast Bio
Form/Appearance	Lyophilized Monoclonal Antibody
Storage	Lyophilized: 2-8°C (desiccated) Reconstituted: Short Term 2-8°C Long Term: -20°C. Avoid repeated freezing and thawing
Note	Reconstitute with Distilled Water. Recommended starting dilution is 1:5,000.
Isotype	IgG1
Clonality	Monoclonal