

# Product Datasheet

## Human TRAIL R3 ELISA Kit 850.840.048

Product Name	Human TRAIL R3 ELISA Kit
Catalog Number	850.840.048
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.31-10 ng/ml
Sensitivity	147 pg/ml
Tested Applications	ELISA
Alternative Names	DcR1, CD263
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity was observed with TRAIL, CD117, IL-6R, IL-2R, CD116, TRAIL R1, TRAIL R2, TRAIL R4, CD178, and Granzyme B.
Specificity	Specific for human TRAIL R3 / DcR1
Sample Types	Serum, Cell lysates