

Product Datasheet

Salmonella paratyphi A LPS Antibody HM514

Product Name	Salmonella paratyphi A LPS Antibody
Catalog Number	HM514
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with other Salmonellae and related gram negative organisms including E. coli and Shigella
Specificity	Salmonella paratyphi A LPS
Matched Pairs	HM039 (capture)
Preservatives	0.1% Sodium Azide