

Product Datasheet

Calicivirus, Feline Antibody HM518

Product Name	Calicivirus, Feline Antibody
Catalog Number	HM518
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	IHC, ELISA
Immunogen	Feline Calicivirus capsid
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG2a, kappa
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	no reactivity with Feline Herpesvirus, FIV or FeLV
Specificity	Feline Calicivirus VP1 capsid protein
Preservatives	0.05% Sodium Azide