

# Product Datasheet

## Human TRAIL R4 ELISA Kit 850.850

|                     |  |
|---------------------|--|
| Product Name        | Human TRAIL R4 ELISA Kit   |
| Catalog Number      | 850.850  |
| Description         | Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range               | 312.5-10000 pg/ml  |
| Sensitivity         | 64 pg/ml   |
| Tested Applications | ELISA  |
| Alternative Names   | DcR2, CD264  |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity was observed with TRAIL, CD117, IL-6R, IL-2R, CD116, TRAIL R1, TRAIL R2, TRAIL R3, CD178 and Granzyme B.   |
| Specificity         | Specific for human TRAIL R4 / DcR2 / TRUNDD  |
| Sample Types        | Serum, Cell lysates  |