

Product Datasheet

Human TRAIL R4 ELISA Kit 850.850

Product Name	Human TRAIL R4 ELISA Kit
Catalog Number	850.850
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	312.5-10000 pg/ml
Sensitivity	64 pg/ml
Tested Applications	ELISA
Alternative Names	DcR2, CD264
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity was observed with TRAIL, CD117, IL-6R, IL-2R, CD116, TRAIL R1, TRAIL R2, TRAIL R3, CD178 and Granzyme B.
Specificity	Specific for human TRAIL R4 / DcR2 / TRUNDD
Sample Types	Serum, Cell lysates