

Product Datasheet

Lincomycin Antibody HM528

Product Name	Lincomycin Antibody
Catalog Number	HM528
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Lincomycin-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, kappa
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Lincomycin 100%; Penicillin, Gentamycin, Tetracycline, Oxytetracycline <1%
Specificity	Lincomycin
Matched Pairs	LA331 (capture)
Preservatives	0.05% Sodium Azide