

Product Datasheet

Salbutamol Antibody HM531

Product Name	Salbutamol Antibody
Catalog Number	HM531
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Salbutamol-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, lambda
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Salbutamol 100%, Clenbuterol and Ractopamine < 1%
Specificity	Salbutamol
Matched Pairs	LA332 (capture)
Preservatives	0.05% Sodium Azide