

Product Datasheet

Dexamethasone Antibody HM533

| | |
|----------------------------|--|
| Product Name | Dexamethasone Antibody |
| Catalog Number | HM533 |
| Description | Ion Exchange Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Immunogen | Dexamethasone-KLH |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1, lamda |
| Clonality | Monoclonal |
| Purity | > 90 % by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl |
| Cross Reactivity | Dexamethasone 100%; DHEA, DHEAs, Cortisol, Corticosterone, Progesterone <1% |
| Specificity | Dexamethasone |
| Matched Pairs | LA333 (capture) |
| Preservatives | 0.05% Sodium Azide |