

Product Datasheet

Dexamethasone Antibody HM536

Product Name	Dexamethasone Antibody
Catalog Number	HM536
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Dexamethasone-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1, lamda
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Cross Reactivity	Dexamethasone 100%; DHEA, DHEAs, Cortisol, Corticosterone, Progesterone <1%
Specificity	Dexamethasone
Matched Pairs	LA333 (capture)
Preservatives	0.05% Sodium Azide