

Product Datasheet

Hepatitis C Virus (HCV-core) Antibody HM540

Product Name	Hepatitis C Virus (HCV-core) Antibody
Catalog Number	HM540
Species/Host	Mouse
Tested Applications	LF, ELISA, WB, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	reactive with genotypes 1a & 2
Specificity	HCV core
Matched Pairs	HM541 (detection)
Preservatives	0.1% Sodium Azide