

Product Datasheet

Hepatitis C Virus (HCV-core) Antibody HM541

Product Name	Hepatitis C Virus (HCV-core) Antibody
Catalog Number	HM541
Species/Host	Mouse
Tested Applications	LF, ELISA, WB, IHC, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	reactive with genotypes 1a, 1b & 2a
Specificity	HCV core (33-38 aa epitope)
Matched Pairs	HM539 or HM540 (capture)
Preservatives	0.1% Sodium Azide