

# Product Datasheet

## Human Perforin ELISA Kit 850.860

Product Name	Human Perforin ELISA Kit
Catalog Number	850.860
Description	Diaclone ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	62.5-2000 pg/ml
Sensitivity	40 pg/ml
Tested Applications	ELISA
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with Trail, TRAIL R3, TRAIL R4, Granzyme B, TNFα, IFNγ, Fas, Fas L, IL-1b, IL-2
Specificity	Recognizes natural human Perforin
Sample Types	Supernatant, Buffered Solutions, Not suitable for use with serum and plasma.