

## Product Datasheet

### Distemper Virus, Canine Antibody HM545

<b>Product Name</b>	Distemper Virus, Canine Antibody
<b>Catalog Number</b>	HM545
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Purified Canine Distemper Virus
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Specificity</b>	Canine Distemper Virus
<b>Matched Pairs</b>	HM547 (detection)
<b>Preservatives</b>	0.09% Sodium Azide