

Product Datasheet

Distemper Virus, Canine Antibody HM547

Product Name	Distemper Virus, Canine Antibody
Catalog Number	HM547
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Canine Distemper Virus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	Canine Distemper Virus
Matched Pairs	HM545 or HM546 (capture)
Preservatives	0.09% Sodium Azide