

Product Datasheet

Pinaca Antibody HM548

Product Name	Pinaca Antibody
Catalog Number	HM548
Description	DEAE Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	DEAE Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	AB-PINACA 5-hydroxypentyl 120%, AB-PINACA 100%, AB-PINACA 5-penanoic acid 80%
Specificity	AB-PINACA
Matched Pairs	LA366 (capture)
Preservatives	0.05% Sodium Azide