

Product Datasheet

Shigella flexneri Antibody HM551

Product Name	Shigella flexneri Antibody
Catalog Number	HM551
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with S. flexneri 1a, 3a, 4a, 5a, 6; Salmonella or E. coli
Specificity	S. flexneri 2a LPS
Preservatives	0.1% Sodium Azide