

## Product Datasheet

### Serum Amyloid A (SAA), Feline Antibody HM558

<b>Product Name</b>	Serum Amyloid A (SAA), Feline Antibody
<b>Catalog Number</b>	HM558
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	WB, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Specificity</b>	Feline SAA
<b>Matched Pairs</b>	HM558 (self-pairs) or HM559
<b>Preservatives</b>	0.1% Sodium Azide