

Product Datasheet

Procalcitonin (PCT) Antibody HM563

Product Name	Procalcitonin (PCT) Antibody
Catalog Number	HM563
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, LF
Immunogen	Recombinant Procalcitonin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	Procalcitonin
Matched Pairs	HM564 (detection)
Preservatives	0.05% ProClin 300