

## Product Datasheet

### Fatty Acid Binding Protein (H-FABP) Antibody HM565

<b>Product Name</b>	Fatty Acid Binding Protein (H-FABP) Antibody
<b>Catalog Number</b>	HM565
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, CLIA
<b>Immunogen</b>	Recombinant Human H-FABP
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	>95%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.0
<b>Specificity</b>	FABP
<b>Matched Pairs</b>	HM566 (detection)
<b>Preservatives</b>	0.02% Sodium Azide