

Product Datasheet

Fatty Acid Binding Protein (H-FABP) Antibody HM565

Product Name	Fatty Acid Binding Protein (H-FABP) Antibody
Catalog Number	HM565
Species/Host	Mouse
Tested Applications	LF, ELISA, CLIA
Immunogen	Recombinant Human H-FABP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Specificity	FABP
Matched Pairs	HM566 (detection)
Preservatives	0.02% Sodium Azide