

Product Datasheet

C-Reactive Protein (CRP) Antibody HM574

| | |
|----------------------------|--|
| Product Name | C-Reactive Protein (CRP) Antibody |
| Catalog Number | HM574 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA, CLIA |
| Immunogen | Purified human CRP |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | > 95 % |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.0 |
| Specificity | Human CRP |
| Matched Pairs | HM573 (capture) |
| Preservatives | 0.02% Sodium Azide |