

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM575

Product Name	Creatine Kinase MB (CK-MB) Antibody
Catalog Number	HM575
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, CMIA
Immunogen	CK-MB recombinant protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 1 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 \pm 0.2
Specificity	Creatine Kinase MB
Matched Pairs	HM576 (detection)
Preservatives	None