

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM576

| | |
|----------------------------|--|
| Product Name | Creatine Kinase MB (CK-MB) Antibody |
| Catalog Number | HM576 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, CLIA, CMIA |
| Immunogen | CK-MB recombinant protein |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term (≤ 1 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 90 % by SDS-PAGE |
| Buffer | 20 mM Phosphate Buffered Saline, pH 7.4 |
| Specificity | Creatine Kinase MB |
| Matched Pairs | HM575 (capture) |
| Preservatives | None |