

Product Datasheet

Chikungunya (CHIKV-Env 2) Antibody HM587

Product Name	Chikungunya (CHIKV-Env 2) Antibody
Catalog Number	HM587
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Chikungunya Virus Envelope 2 Protein
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a, kappa
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Dengue, Zika, YFV, WNV, TBV, JEV
Specificity	CHIKV E2
Preservatives	0.05% Sodium Azide