

Product Datasheet

Human IL-10 HS ELISA Kit 850.880.192

| | |
|---------------------|---|
| Product Name | Human IL-10 HS ELISA Kit |
| Catalog Number | 850.880.192 |
| Description | Diaclone High Sensitivity ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure. |
| Range | 1.56-50 pg/ml |
| Sensitivity | 0.98 pg/ml |
| Tested Applications | ELISA |
| Alternative Names | NULL |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human cytokines |
| Specificity | Recognizes both natural and recombinant human IL-10 |
| Sample Types | Serum, Plasma, Cell culture supernatant |