

Product Datasheet

Pinaca Antibody HM594

Product Name	Pinaca Antibody
Catalog Number	HM594
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Pinaca-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	AB-PINACA, AB-PINACA 4-hydroxypentyl & AB-PINACA 5-hydroxypentyl 100%; AB-PINACA 5-pentanoic acid 95%; APINACA & APINACA 5-hydroxypentyl 1%
Specificity	AB-PINACA
Matched Pairs	LA384 (capture)
Preservatives	0.05% Sodium Azide