

Product Datasheet

UR-144 Antibody HM595

Product Name	UR-144 Antibody
Catalog Number	HM595
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	UR-144-BSA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	UR-144 100%; UR-144 5-hydroxypentyl, UR-144 5-pentanoic acid, UR-144 4-hydroxypentyl = 110%
Specificity	UR-144 and UR-144 metabolites
Matched Pairs	LA385 (capture)
Preservatives	0.05% Sodium Azide