

Product Datasheet

Pinaca Antibody HM616

Product Name	Pinaca Antibody
Catalog Number	HM616
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	APINACA 5-hydroxypentyl 100%; APINACA 80%; UR-144 4-hydroxypentyl 10%; UR-144 5-hydroxypentyl; UR-144, UR-144 5-pentanoic acid, AB-PINACA, AB-PINACA 5-hydroxypentyl & AB-PINACA 5-pentanoic acid = 3%
Specificity	APINACA
Matched Pairs	LA366 (capture)
Preservatives	0.05% Sodium Azide