

# Product Datasheet

## Human IFN gamma HS ELISA Kit 850.900.048

Product Name	Human IFN gamma HS ELISA Kit
Catalog Number	850.900.048
Description	Diaclone High Sensitivity ELISA Kits include pre-coated strip plates, biotinylated secondary antibody, standards, controls (where applicable), buffers, streptavidin-HRP, TMB, stop reagents, and a detailed procedure.
Range	0.78-25 pg/ml
Sensitivity	0.69 pg/ml
Tested Applications	ELISA
Alternative Names	IFN-g,IFN-g, IFN-gamma IFN-gamma
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IFN-g
Sample Types	Serum, Plasma, Cell culture supernatant