

## Product Datasheet

### Anti-Mullerian Hormone (AMH) Antibody HM621

<b>Product Name</b>	Anti-Mullerian Hormone (AMH) Antibody
<b>Catalog Number</b>	HM621
<b>Description</b>	Protein A Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	LF, ELISA, CLIA
<b>Immunogen</b>	Recombinant AMH
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Note</b>	Limit of Detection for Lateral Flow 1 ng/ml (colloidal gold) and ≤ 20 ng/ml fluorescent label
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Specificity</b>	AMH N-terminal domain
<b>Matched Pairs</b>	HM620 (capture)
<b>Preservatives</b>	0.05% ProClin 300