

Product Datasheet

Anti-Mullerian Hormone (AMH) Antibody HM621

Product Name	Anti-Mullerian Hormone (AMH) Antibody
Catalog Number	HM621
Description	Protein A Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA, LF
Immunogen	Recombinant AMH
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Limit of Detection for Lateral Flow 1 ng/ml (colloidal gold) and ≤ 20 ng/ml fluorescent label
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	AMH N-terminal domain
Matched Pairs	HM620 (capture)
Preservatives	0.05% ProClin 300